


# LUNG BIOMARKERS FOR TOXICOLOGY STUDIES



**Joe Mauderly**

**Lovelace Respiratory Research Institute  
Albuquerque, NM**



# “BIOMARKERS” DISCUSSED

Respiratory function

Lung imaging

Bronchoalveolar lavage

Tissue chemistry

Special stains & immunohistochemistry

Cell proliferation

Morphometrics

Gene expression

Definition

Uses

Advantages

Disadvantages

## Suitability for complementing NTP hazard assessment bioassays



- Useful for detecting and characterizing adverse responses
- Acceptably standardized methods
- Generally accepted interpretation of results & extension to humans
- Enhances evaluation of hazard by current definitions

# RESPIRATORY FUNCTION


## Ventilation

Respiratory frequency =  $f$   
 Tidal volume =  $V_T$   
 Minute volume =  $V_E$


## Mechanics



## Volumes




## Diffusion and Perfusion



# RESPIRATORY FUNCTION

## Airway responsiveness



Singh et al., *Am. J. Resp. Crit. Care Med.* 168: 342, 2003



# RESPIRATORY FUNCTION

## Uses

- Describes functional manifestation of structural changes
- Places functional impacts of disorders into clinical context
- Provides correlates to humans

## Advantages



- Can be done non-destructively (with anesthesia)
- Established methods and interpretation

## Disadvantages

- Most tests require specialized equipment and expertise
- Not very sensitive to subtle or scattered tissue changes
- Not a substitute for histopathology or other indicators
- Does not determine type of morphological abnormality

# LUNG IMAGING

Gas density images of rat lungs using MRI measurement of apparent gas diffusion coefficient of hyperpolarized He



Diffusion coefficients of middle third of lungs

Chen et al, Proc. Nat. Acad. Sci. 97: 11478, 2000

Image lungs using x-ray, MRI, PET scan, etc.

With or without contrast media to highlight structures

# LUNG IMAGING

## Uses

Detection and staging size and distribution of lesions

Following abnormalities with time

Provides correlates to humans

## Advantages

Non-destructive (requires anesthesia)

## Disadvantages


Marginal resolution of some methods for imaging rodents (e.g., x-ray)

Some methods require sophisticated/expensive equipment (e.g., MRI)


Not very sensitive to subtle or scattered, focal tissue changes

# BRONCHOALVEOLAR LAVAGE

Biomarkers in fluid instilled into and withdrawn from airways



Henderson, *Exper. Toxicol. Pathol.* 57: 155, 2005



Mauderly, *Lab An. Sci.* 27:255, 1977

**Cells** (RBC, WBC, differentials, collections for cell morphology, DNA, or function)

**Protein** (total, albumin, hemoglobin, proteomics)

**Enzymes** (Lactate dehydrogenase,  $\beta$ -glucuronidase, alkaline phosphatase)

**Cytokines/chemokines** (interleukins,  $\text{TNF}\alpha$ , MIP-2, etc.)


**Growth factors**

**Glutathione** (total, reduced)

**Fibronectin**

**Antibodies**

**Elastin/collagen breakdown products**



# BRONCHOALVEOLAR LAVAGE

## Uses

Detection of inflammation, cytotoxicity, oxidative stress, lung tissue metabolism, allergic reactions

Provides correlates to humans

## Advantages

Readily done at necropsy (can be done on one lobe)

Assays are straightforward (kits available for many variables)

Established methods and interpretation

Non-destructive to tissue (can be coupled with other assays)

## Disadvantages


Reflects responses in airspaces, not necessarily in tissues

Modest ability to detect or stage chronic lung disease (e.g., fibrosis, emphysema, cancer)

All species can be lavaged in vivo, but not all survive or completely clear atelectatic spots (mice, S. hamsters, rabbits, NHP, dogs do – rats, guinea pigs, & gerbils do not)


# TISSUE CHEMISTRY

**Lung Collagen  
Radiation Pneumonitis-Fibrosis**



Pickrell et al., *Rad. Res.* 74:363, 1978

**HO-1 in Lung  
With and Without Diesel Emissions Control**



McDonald et al., *EHP*, 112: 1307, 2004

**Chemical analysis to determine content of target material**

**Structural (e.g., collagen) or functional (e.g., hemoxygenase-1, glutathione) analytes**

# TISSUE CHEMISTRY

## Uses

- Measure changes in concentration of specific structural components

- Measure functional status

## Advantages

- Provides quantitative information

## Disadvantages


- Isolation to particular anatomic site is dependent on level of tissue dissection (whole lung, lobe, microdissection)

- Averages all structures and cells within sample


- Destructive

# SPECIAL STAINS AND IMMUNOHISTOCHEMISTRY

## Alcian Blue/Periodic Acid Schiff Stain




**Mucosubstance Volume in Airways of Mice Exposed to Cigarette Smoke**



March *et al.*, *Toxicol. Sci.* 92:545, 2006

## Phosphorylated IKK $\alpha$ Stain (inhibitor of nuclear factor K $\beta$ kinase)



**Activation of NFK $\beta$  in airway of mouse infected with cowpox virus**

J. Hutt, LRRI, personal communication, 9/15/06

**Stain bound chemically (collagen, elastin, mucus, etc.)**

**Stain linked to antibody (BrdU, Bcl-2, cytokeratin, NFK $\beta$ , etc.)**



# **SPECIAL STAINS AND IMMUNOHISTOCHEMISTRY**

## **Uses**

Enhance identification, visualization, and quantitation of structural/chemical abnormalities

## **Advantages**

Provides information on anatomic location and variation

Can target very specific materials

Can be quantitative

## **Disadvantages**


Reliant on sampling/sectioning scheme

Provides 2-dimensional view

Immunohistochemistry depends on availability of antibody

# CELL PROLIFERATION

**BrdU staining (brown) of airway epithelial cells**



Y. Tesfaigzi, LRRI, personal communication, 9/15/06

**Identify dividing cells using nuclear marker**

**Can label ante-mortem (e.g., BrdU - thymidine analog) or post-mortem (e.g., Ki67 - nuclear antigen)**

**Often measured as cells/mm basal lamina**

# CELL PROLIFERATION

## Uses

- Assess cell turnover rates

- Detect and measure proliferative events

## Advantages

- Provides quantitative information on proliferative status

- Detect “harbinger” cellular responses

## Disadvantages

- Ante-mortem labeling by injection provides limited time window

- Prolonged ante-mortem labeling requires repeated injections or implantation of osmotic pumps


- Concurrent markers often necessary to confirm cell type

- Manual counting takes time


- Proliferation rate may be only slightly above background in chronic conditions (may miss time course of early events)

# MORPHOMETRICS

Horizontal Cycloid Grid Overlying Alveoli



Intercepts with tissue



Width of alveoli

## Air space enlargement

$$L_m = \text{length} \div \sum I_{\text{spt}}$$

$$V_{\text{Vair}} = \sum P_{\text{air}} \div \sum P_L$$

## Tissue destruction

$$V_{\text{Vspt}} = \sum P_{\text{spt}} \div \sum P_L$$

$$S_a = (4 \times V_L) \div L_m$$

T. March, LRRI, personal communication, 9/15/06

## Quantitative analysis of structural dimensions

- 1) Extent and distribution of abnormalities (e.g., collagen)
- 2) Morphological dimensions (e.g., airspaces, airway and vessel walls)

**Stereology** (extrapolation to 2-dimenional surface density or 3-dimensional volume density by 2-dimensional point counting)

**Linear frequency rates** (counts of events per unit distance, e.g., cells per mm basal lamina)

# MORPHOMETRICS

## Uses

Provides quantitative assessment of structural changes

## Advantages

Adjunct to subjective, non-quantitative assessment/scoring

Supports statistical analyses

## Disadvantages

Can be tedious

Accuracy highly dependent on:

1) identification of structures


2) sampling and computational schemes

number & location of sections

number of events counted

accuracy of normalization variable (e.g., lung volume)

# GENE EXPRESSION



Liekauf, presented at NERC Annual Meeting, May 2003, unpublished

**Expression of genes, using mRNA isolated from lung**

**rtPCR/gel or real-time rtPCR for specific genes**

**Micro-array survey (standard or custom arrays)**

**Comparison with control expression level**

# GENE EXPRESSION

## Uses

- Assess activation status of specific gene(s)
- Screen for detecting and categorizing exposure effects
- Generate hypotheses about mechanisms of response

## Advantages

- Sensitive – likely to detect exposure-related differences
- Patterns may be descriptive of response type and mechanisms

## Disadvantages

- Usually requires fresh or frozen tissue
- Integrates across all cell types in sample
- Gives information for one point in time
- Gene expression may not reflect protein production
- Expression of “clock” genes is affected by sample time
- Micro-arrays incur data analysis challenges
- Utility dependent on understanding links between gene expression and outcomes relevant to human hazard (limited number of well-established extrapolation pathways)

# PROTEOMICS

## Uses, Advantages, Disadvantages

Generally same as for gene expression

Advantage: measures actual product

Advantage or disadvantage: measures total product present,  
regardless of when produced



# SUMMARY

- **All of these biomarkers (and others) could be useful**
- **Probably none should be included routinely in all bioassays**
- **These are not substitutes for conventional bioassay endpoints**
- **Most of these are not high-throughput approaches**
- **Selection depends on:**
  - 1) Outcome of concern
  - 2) Nature of evidence for human hazard
  - 3) Current regulatory definitions of “hazard” and “adverse effect”
  - 4) Extent to which mechanisms and time course are of interest