

---

# Concept Review:

---


## Analytical Chemistry for the Environmental Toxicology Program



# Why is Analytical Chemistry Important in Bioassays?

---

## Bioassay Paradigm



Response must be supported with good analytical chemistry data for exposure and internal dose.

# Classes of Test Articles Studied by the NTP

---

Drinking Water Contaminants - BDCM, DCA, Bromopicrin

Food Additives and contaminants — Methyleugenol, hydroxymethylfurfural, 2- and 4-Methylimidazole

Materials used in making plastics - bis(2-Chloroethoxymethane), Dimethyl-p-toluidine, Formamide

Flame Retardants — TBBPA, PBDEs

Consumer Product Ingredients — Dibromodicyanobutane, Chitosan

Metals - Hexavalent chromium, Chromium picolinate

Pharmaceuticals - AIDS Combination Therapies, Elmiron

Solvents - MIBK, Tetralin, Stoddard Solvent

Botanical Products - Ginkgo biloba extract, Goldenseal root

Blasting Agents - Blasting sand, Garnet, Crushed glass, Coal slag

# Overview of the Chemistry Role

---

## Types of studies supported:

- Carcinogenicity
- General Toxicology
- Reproductive Toxicology
- Immunotoxicology
- Genetic Toxicology
- Differential Gene Expression
- DIR in-house research upon request

## Overview of Chemistry Role (cont'd)

---

### Types of tasks performed:

- Chemical Procurement
- Chemical Characterization
- Dose Formulation Development
- Biological Sample Analysis
- Toxicokinetics Studies with Unlabeled Compounds

# Capabilities

---

## Chemical Characterization:

- Physical Constants Determination
- IR, NMR, MS, MS/MS
- Elemental Analysis
- Water Determination
- Chromatographic Analyses - HPLC, GC, IC, TLC – all detectors
- Functional Group Titration
- Storage Stability Evaluation

# Chemical Characterization for Bioassays

---

## Unequivocal Identity

- IR, NMR, MS
- Physical Constants (Chronic only)

## Purity Determination

- Water Determination
- Elemental Analysis (Chronic only)
- 2 Orthogonal Chromatographic Analyses (Organics, ICP/AES or ICP/MS for Inorganics)
- Impurity identifications at  $\geq 1$  %
- Impurities reported at  $\geq 0.1$  %

# Preliminary Chemistry Studies

---

- Solubility
- Suspendability
- Palatability
- Gavageability
- Inhalation Feasibility



# Dose Formulation

---

- Vehicle Characterization
- Dose Analysis Method Development and Validation
- Homogeneity Evaluation
- Stability Studies

# Toxicokinetics Studies

---

- Unlabeled test article
- Pilot study – feasibility, when needed
- Preliminary study - IV and bioassay route
- Initial doses based on literature values of LD<sub>50</sub>
- Tissues used to develop bioanalytical method
- Definitive study - GLP with multiple samples/time point
- Single animal data QC'd spreadsheets => modelers for PBPK
- Results => parameter values from non-compartmental model

# Typical Non-compartmental Parameters

---

Values are reported for typical parameters  
including:

- AUC
- AUMC
- CI
- F
- MRT
- t
- $\bar{V}$
- k

# New to this SOW

---

Current impurities paradigm = Identify >1 %; report >0.1 %

Some studies require more:


## New Assignment – Low Level Impurity Determination (LLID)

- Offline from routine characterization
- Iterative plan with interim data submitted
- High priority

**Example: Anthraquinone**

# Anthraquinone Purity Analysis by GC/FID

---




Initial analysis by GC/FID - 0.1 % impurity

Others - smaller and therefore not reported

Overall purity = 99.9 %

# Anthraquinone Purity Analysis by HPLC/UV

---



Initial analysis by HPLC/UV – reported impurities were 0.2 and 0.3 % of total peak area

Others - 0.01 % and 0.03 % not reported

Overall purity = 99.5 %

# Impurity Questions

---


In response to concerns of stakeholders:

- 1) Impurities identified
- 2) Rationale - difference between GC purity and HPLC purity
- 3) Impurities quantitated


# Impurity Questions 1) and 2)

---


1) HPLC/MS showed the presence of 5 compounds:




Anthraquinone




Anthracene



9-Nitroanthracene



Phenanthrene




Anthrone

2) Ultraviolet absorbance roughly doubles with each conjugated double bond



# Impurity Question 3)

## Anthracene Impurity Quantitation by Standard Addition with GC/FID



Anthracene = 0.05 %


9-Nitroanthracene = 0.1 %

Phenanthrene, Anthrone = < 0.002 %

Overall purity = 99.85 %

## Impurity Question 3) cont'd

Quantitation of Anthracene Impurity with HPLC/UV  
By Standard Addition



Anthracene = 0.06 %

9-Nitroanthracene = 0.11 %

Phenanthrene, Anthrone = < 0.001%

Overall purity = 99.83 %

# New to this SOW

---

## New assignment – Chemical Identity and Purity Screen (CIPS)

Acceptable techniques will be:

- Automated
- Easily interpreted
- Cost effective

Chief techniques will be NMR (organics) and ICP/AES (inorganics). Fall-backs will be Flow-injection or Direct Probe MS and HPTLC

# New to this SOW

---

## New Assignment – Biochemical Measurement (BCM)

Intended as an interim measure.

Assignments will be:

- documented in the literature
- well accepted in the research community
- require no method development
- biochemical measurements

Examples – chemical biomarker measurements, enzyme assays, protein binding assays

# Expectations for Procurement

---

- 3 awards projected
- 800 assignments/year/award
- Full range of capabilities/award

# Cost Containment

---

1. Use test articles/methods in multiple programs wherever possible
2. Use high throughput analyses for test articles intended for high throughput studies
3. Place emphasis on reporting interim data to direct work in a facile timeframe