

Nomination of the whole blood pyrogen test for extension to non-endotoxin pyrogens

Thomas Hartung

Doerenkamp-Zbinden Professor and Chair for Evidence-based Toxicology,

Director, Center for Alternatives to Animal Testing (CAAT)

Joint appointment: Molecular Microbiology and Immunology

Johns Hopkins University, Baltimore, US


Professor of Pharmacology and Toxicology, University of Konstanz, Germany



Pyrogen (bacterial and fungal contamination) testing still consumes more than 300.000 rabbits per year



Limulus assay is restricted to Gram-negative endotoxin, disturbed by many substances and does not reflect human pyrogen potency



**1995 Whole blood pyrogen test
2003 and 2004 Validated
2006 Validity statement
2009 Accepted Eur. Pharmacopoeia
and with limitations by ICCVAM / FDA
PyroDetect (Biotest)**

Novel Pyrogen tests based on the human fever reaction

Exogenous
Pyrogen



1. Whole Blood
2. Leukocytes
3. Cell lines
4. Cryopreserved blood cells



Endogenous
Pyrogen




Fever

ELISA

The rabbit test as reference

The temperature reaction of a sensitive strain (n = 171)

- limit temperature increase in Europe: 0.55° C
 ➔ 5 EU/10 ml = 0.5 EU/ml defined as threshold
- 13 drugs with five contaminations each (0, 0.25, 0.5, 1 EU/ml)



challenging original study design

➔ (58% specific, 88% sensitive for rabbit)

Results – Predictive Capacity

Joint Research Centre

Prediction Model 'PM0'

- non-pyrogenic: 0 EU/ml
 - pyrogenic: 0.25, 0.5, 1 EU/ml
- ➡ 30 non pyrogenic and 120 pyrogenic samples

Original Prediction Model

- non-pyrogenic: 0, 0.25 EU/ml
 - pyrogenic: 0.5, 1 EU/ml
- ➡ 60 non pyrogenic and 90 pyrogenic samples

Test	Specificity	Sensitivity	Test	Specificity	Sensitivity
PBMC	96.7%	100%	PBMC	93.1%	77.6%
WBT A	92.6%	96.2%	WBT A	98.8%	83.6%
WBT B	85.2%	99.0%	WBT B	82.4%	89.1%
WBT fresh	85.7%	99.1%	WBT fresh	97.4%	81.4%



Conclusion

The validated methods using cryopreserved cells

- are an important standardisation
 - ➔ widely applicable by overcoming limitations
 - availability of primary cells
 - exclusion of abnormal donors
 - test for infectious agents
- are reproducible and robust
- perform well in terms of predictive capacity
- should augment regulatory acceptance

***Adverse fever reactions by an
infusion solution containing gelatine
negative Limulus Test as release criterion
problematic batches recalled***

samples blinded and sent to PEI

fever	rabbit	limulus	IL-1 pg/ml	IL-6 pg/ml	TNF pg/ml
-	-	-	8.5	28.0	28.2
+	+	-	142.6	654.4	67.6
+	-	-	421.5	9444.0	116.7



Non-endotoxin stimuli in the different tests

	pyrogen	LAL	blood-IL1	blood-IL6	PBMC-IL6	MM6-IL6	THP-1-Neo
Curdian	-	+	-	-	-	-	-
Glucan-Barley	-	+	-	-	-	-	-
Glucan-Yeast	-	+	-	-	-	-	-
Glucan STD	-	+	-	-	-	-	-
Lipid A	-	+	-	-	-	-	-
PHA-E	-	+	-	-	-	-	-
PHA-L	+	+	+	+	+	+	+
Zymosan	-	+	-	-	-	-	-
LPS	+	+	+	+	+	+	+
LTA	+	-	+	+	+	+	+



New challenges

Patient safety: non-endotoxin pyrogens

**New expression systems in gene technology
(Gram-positive bacteria, insect and fungal cells)**

Medical devices

Cellular therapies


Environmental pyrogens

ICCVAM Recommendations: Future Studies May 2008

1. Both endotoxin-spiked and non-endotoxin spiked samples should be included. Non-endotoxin standards should be characterized prior to their use in any study, if possible.
2. ... Good Laboratory Practice.
3. ... products that have intrinsic pro-inflammatory properties ...
4. Optimally, a study that includes 3-way parallel testing, with the *in vitro* assays being compared to the RPT and the BET, should be conducted to allow for a comprehensive evaluation of the relevance and comparative performance of these test methods. These studies may be conducted with historical RPT data provided that the same substances (i.e., same lot) are tested in each method. Based on ethical and scientific rationale, any *in vivo* testing should be limited to those studies that will fill existing data gaps.
5. Test substances that better represent all categories of sample types (e.g., pharmaceuticals, biologicals, and medical devices) intended for testing by the methods should be included.

Bacterial pyrogens

Gram-negative
bacterial
cell wall




LPS

(Lipopoly-
saccharide)

Medline publ.
~ 50 000

Gram-positive
bacterial
cell wall



LTA

(Lipoteichoic
acid)

Medline publ.
~ 1 000

LTA is a Gram-positive endotoxin

Contra

inactive when purified

LPS contamination

synthetic LTA inactive

less potent than LPS

No specific inhibitor

Pro

similar structure

conserved principle

immunostimulatory

LTA is a Gram-positive endotoxin

Contra

~~inactive when purified~~

Morath 2001, J. Exp. Med.

~~LPS contamination~~

Morath 2001, Inf. Immun.

~~synthetic LTA inactive~~

Morath 2002, J. Exp. Med.

Deininger 2003, J. Immunol.

Deininger 2007, Clin. Vacc. Imm.

von Aulock 2003, Immunobiol


Deininger 2008, Immunobiol

Draing 2008, Eur. J. Immunol.

~~less potent than LPS~~


~~No specific inhibitor~~

Comparison of different donors




➡ Different absolute but similar relative response of donors

Relative response of donors



Differences caused mainly by different monocyte counts

Surface bound LTA is more potent



➔ Presentation of LTA on a surface increases 1000fold its inflammatory potency

Deininger, Immunobiol. 2008


Systematic review following Evidence-based Medicine as to the fulfillment of Koch/Dale criteria for the major candidate structures inducing cytokine release in human monocyte / macrophages

Koch / Dale	LTA	PGN	Lipoproteins
Presence in bacteria	+++	ubiquitous	-
Synthesis inhibition abrogates cytokines	+	not possible	+
Isolated substance induces cytokines	+++	++	++
Block of substance	+	-	-

1-5 publications: +; 6-10 publications: ++; >10 publications: +++

Rockel, 2011, submitted

Why is the nature of the Gram-positive endotoxin 50 years after identification of the Gram-negative one still under discussion?



- **LTA instability**
- **Presentation effect (potency)**
- **Differences mouse to human**
- **Not complete endotoxin (IFN γ)**
- **Limited commercial availability**
- **No specific test comparable to LAL**

Toward a follow-up validation

- **Set up a validation management group**
- **Design the validation study**
- **Choice of test materials and non-endotoxin pyrogens**
- **GLP?**
- **Suggested:**
 - **LTA and lysates of Gram-pos bacteria**
 - **Rabbit testing only for spike materials**
 - **Limited reproducibility re-assessment**